

## Artificial intelligence augmented ultrasound detection of hip dysplasia

Jacob Jaremko, MD, PhD  
MALMIC April 2023

1

## My current roles

- Professor, Radiology
- Adjunct Professor, Computing Science
- Canada CIFAR AI Chair
- Partner, Medical Imaging Consultants
- Co-founder, MEDO.ai
- Fellow, AMII

2

## Developmental Dysplasia of Hip (DDH)

- Common (1/100)
- Abnormal bony shape → unstable joint → premature arthritis
- Preventable if detected in infancy

3

## Original idea Down Under

- During clinical fellowship training
- Learning to perform hip dysplasia ultrasound
- Sonographer Cain Brockley: why can't we look at whole 3D shape of hip?

4

## DDH Screening

- Every day in Alberta...
  - ~4 children (1 Indigenous) born with DDH
  - ~2 hip arthroplasties for OA from DDH
- 95% of infant DDH is cured by wearing a soft harness for 6 weeks
- Graf ultrasound screening
  - Many false positives
  - Too costly for general use
- Can we screen with US/AI?

5

## DDH – Alpha Angle

6

## Ideal scenario for screening

- Relatively common disease
  - Especially in Indigenous
- Costly lifelong morbidity when missed
- Highly effective treatment available when detected.
- Harmless screening test (ultrasound)
- Why don't we screen worldwide?

7

## 2D US not recommended for screening

- False positives
  - Recall for costly follow-up exams
- False negatives
  - Missed late presentation DDH
- Unreliable indices
- Limited prognostic information
  - Cochrane review 2011, AAP (Shaw 2016), "Choosing Wisely"

The Opportunity:  
Screen for DDH with a better test!

8

## Hip 3D Ultrasound

Commercial linear probes  
Mechanical sweep  
Scan time:



Philips 3.2 s, Xmatrix 0.6 s  
Canon <2 s



9

## 'Sweep' video ultrasound

- User saves the entire video clip as they 'sweep' through the hip.
- Possible on any modern ultrasound probe
- Captures the whole hip shape
- Expert user or AI can analyze images afterward: normal or dysplastic?
- "3D" ultrasound = automated sweep





10

**Radiology**  
Doroosh Zoroobi, MEng, PhD  
Aljilash Harondranathan, PhD  
Esmail Mostafaei, BSc  
Myles Mabey, MSc  
Saba Patha, PhD  
Diana Cobzas, PhD  
Padma Rao, MSc, MRCP, FRCP, FRACR, FRANZCR  
Sukdeep K. Duttal, MD, BSc, MRSc, FRCS, FRCS  
Jeevesh Kapur, MD, MMed, FRCP  
Jacob L. Jaramico, MD, PhD, FRCP

**Developmental Hip Dysplasia  
Diagnosis at Three-dimensional  
US: A Multicenter Study<sup>1</sup>**  
Radiology: Volume 287: Number 3—June 2018 • radiology.rsn.org


- We tested 'sweep' method using 3D ultrasound probes
- Borderline = Graf IIa (return for follow up scan)
- In 1,728 hips, 3DUS reduced **borderline** cases by nearly **75%**



11

## 3D ultrasound image analysis


- Manual segmentation
- Tedious - 250 slices
- "Labor of love"
- Can we automate this?



12

## Easy screening?


- If analysis could be automated → population screening



13

## Artificial Intelligence (AI) Image Analysis


- Smartphones now recognize faces
- Computers can now also recognize pathology in medical images.
- Convolutional neural network (CNN)
  - Input: entire image
  - Output: *numbers, or an image highlighting abnormalities.*



14

## AI software to assist radiologists

- Example: Osteodetect
- US FDA approved: finds fractures in Xray by AI
- Many similar products rapidly emerging



15

## EU Approves First AI For Analyzing X-Rays Automatically

By Adrianna Nine on April 7, 2022 at 5:30 am [Comments](#)

- Oxipit ChestLink; has CE mark
- 3-step onboarding per site:
  - Initial retrospective phase
  - Prospective 'in parallel' phase
  - Autonomous phase
- All 'normal with high confidence' CXR can be removed from workflow (expected 15-40% of CXR)
- Heatmaps supplied for abnormal CXR



16

## ChestLink performance

- Tested in Finland
- 10,000 CXR, clinics; 45% normal
- AI picked 1,692 normal, of which 9 false negative
  - ?cardiomegaly, ?atelectasis, 1 small effusion
  - No serious miss; SN 99.8%, SP 36.4%
- AI could have removed 1/3 of normal CXR from queue
- Cost savings? Relevance to hospital work? Benefits?
  - Keski-Filppula 2022, unpublished




False negative: ?L effusion


17

## AI Strategies

- CNN → Shape analysis
  - identifies anatomy; output is shape mask
  - Then conventional methods for indices & diagnosis




- End to end ANN
  - Initial layers are convolutional
  - Final fully connected layers predict Dx



18

## Automated image segmentation

- Deep learning to replace user identification of seed points
- Series of 2D and 3D CNN, use info from adjacent slices
- Eventually successful



19

## Interest is spreading

- Vancouver



- China



20

## Automated image segmentation

- Getting easier on natural images
- 2023: Facebook Segment-Anything model
- Likely also works with ultrasound



21

## MEDO.ai



- Med-tech startup company
- Singapore / Edmonton 2018
- Acquired by Exo Inc. (USA) in mid 2022



22

## MEDO Web App (2019)

- Works on browser, eg. Chrome
- Automatically detects hip anatomy
- Calculates probability of hip dysplasia
- Generates report to assist radiologist
- US FDA approved; Health Canada ITA granted.



23

## Tablet/Smartphone App (2021)

- US probe plugs in via USB
- App runs entirely on an Android OS device
- No delays for upload/download
- Near-instant diagnosis
- Can use at any point of care



24

## Tablet/Smartphone App (2021)

- Video of tablet screen while scanning baby




25

## Mobile app report

26

## DDH Screening Pilot Study for AICE


- 3 Sites
  - Spruce Grove
    - Westland
    - Westview
  - Red Deer
    - Saint Mary's
- Supported by grant (AI-CE)



27

## Results

- 2 years (2021-23)
- 799 scans / 679 infants
- 90% scanned only by nurses
- 11 babies (1.5%) DDH
  - No risk factors or signs: 5
  - Indigenous: 1
  - 3 cast, 2 surgical
- Internal FU in 52 (8%)
  - 26 AI+, 26 technical failure
- External FU in 23 (3%)
  - SP 98%, PPV 48%
  - SN, NPV TBD...



	+ DDH	- DDH	
+ Scan	11	12	PPV = 48%
- Scan	0	674	-
			Specificity = 98%

28


## Pilot study conclusions

- US/AI DDH screening at point of care is feasible
- Results ~ expert screening
  - Rosendahl 1994 in-hospital: 14% FU rate, ours 3-8%
  - DDH rate detected = 1.5%
- Already helping babies
  - 5 treated who would have been missed
- Next steps:
  - Integrate hardware into app
  - Add treatment at point-of-care
  - Implement zone-wide, Indigenous sites

29

## AI in ultrasound - opportunity

- Handheld ultrasound at point of care
- Don't bring patient into hospital, bring expertise to the patient
- AI to interpret images for users
- 21<sup>st</sup> century stethoscope
- Example of successful commercialization in progress



30

